

EASY GRIP DIGITAL SOLUTIONS

IMPLANT LINE	DIGITAL LIBRARY			PROSTHETIC COMPONENTS					SOFTWARE CAD-CAM
	FOLDER	NAME	PROSTHETIC SOLUTIONS	SCANBODY			ABUTMENT	DIGITAL REPLICA	
				CODE	DESCRIPTION	LASER			
Easy Grip HEX Series 19 - 29 - 39 - 59 - 69	TFI System EasyGrip HEX 7.0	Easy Grip HEX Bridge	Bar prosthetic / Toronto	SBFX Titanium	Ø4.25 2 planes H 5.0	-	PTS00	LD3	EXOCAD 3SHAPE DENTAL WINGS HYPDOCAD
		Easy Grip HEX TBase	Prosthesis on indexed T-Base				PTS34, PTS34-2/4/6		
			Prosthesis on non-indexed T-Base				PTR34, PTR34-2/4/6		
		CONNESSIONE	Indexed prosthesis with direct connection				-		
		CONNESSIONE ROTANTE	Non-indexed prosthesis with direct connection	-					
	TFI System EasyGrip HEX MUA 7.0	CONNESSIONE	Direct screwed prosthesis on straight and angled MUA	SBPS Peek or Titanium	Ø4.50 2 planes H 7.0	-	PS34-0/1/2/3/4/5	LDS4	
		Link4	Screwed prosthesis on straight MUA with link				PS322/24/30/34	LDS5	
		Link5	Screwed prosthesis on angled MUA with link				LPS4	LDS4	
			LPS5				LDS5		
Easy Grip CONE Series 79	TFI System Easy Grip CONE TBASE	LINK	Prosthesis on indexed T-Base	7SBFX Titanium	Ø4.25 2 planes H 5.0	T	7PTS41/2/3	7LD	EXOCAD 3SHAPE DENTAL WINGS HYPDOCAD
		LINK ROT	Prosthesis on non-indexed T-Base				7PTR41/2/3		
		CONNESSIONE	Indexed prosthesis with direct connection				-		
		CONNESSIONE ROTANTE	Non-indexed prosthesis with direct connection				-		
Easy Grip CONE e NARROW Series 79 e 89	TFI System Easy Grip CONE MUA	CONNESSIONE	Direct screwed prosthesis on straight and angled MUA	7SBPS Titanium	Ø4.80 2 planes H 5.0	M	7PS4-1/2/3/4 7PS17-3/5 7PS30-3/5	7LDS	EXOCAD 3SHAPE DENTAL WINGS HYPDOCAD
		LINK	Screwed prosthesis on MUA with indexed link				7LPS		
		LINK ROT	Screwed prosthesis on MUA with non-indexed link						
Easy Grip NARROW Series 89	TFI System Easy Grip NARROW	LINK	Prosthesis on indexed T-Base	SBSX Titanium	Ø4.00 2 planes H 5.0	-	PTS24-2/4/6	LDM	EXOCAD 3SHAPE DENTAL WINGS HYPDOCAD
		LINK ROT	Prosthesis on non-indexed T-Base				PTR24-2/4/6		